

## Transformer Scooter

- No Tools Required
- Folds Compact & Portable
- Up to 35 Miles Per Charge
- Lightweight
- Easy to Carry
- 300lb Weight Capacity
- Remote Operated
- Ideal for Travel
- Fits in Most Vehicles
- Airline Approved Battery



Shown in Yellow  
Available in Four Colors



Specifications	Transformer Scooter
Overall Length	37.4"
Overall Width	22.25" w/ Armrests
Tires	Front – 6x1.5" Rear – 7x2.36"
Maximum Speed	3.7mph
Range	Up to 13.5 miles
Turning Circle	55"
Weight (w/ battery)	53lbs
Color	Red, Blue, Yellow, Black
Brake System	Intelligent Regenerative Electromagnetic Braking
Drive System	Rear Wheel Drive
Handlebar	PU Foam
Weight Capacity	300lbs
Seat Width	17"
Seat Height (from ground)	21.75"
Battery	24V10Ah Lithium
Motor	120W
Battery Charger	DC24V/2A
Folded Dimensions	17" x 22.3" x 25"



1/2021

