



# TRIAXE >>>

## Sport

- Aluminum Frame
- Folds Flat
- Front Wheel Drive
- 44lbs (without battery)
- Dual Rear Disc Brakes
- Speeds to 12mph
- 350lb Weight Capacity
- Battery - 48V 12amp (48V 6amp optional)
- Twist grip Acceleration
- Up to 35 miles per Charge
- 34" long by 20 1/2" wide
- 13" seat (optional 18")
- 6 Limited Edition Colors

## Cruze



- Aluminum Frame
- Folds Flat
- Rear Wheel Drive
- 58lbs (without battery)
- Electromagnetic Braking
- Speeds to 6mph
- 350lb Weight Capacity
- Battery - 24V 12ah
- Paddle acceleration
- Up to 18 miles per Charge
- 39" long by 19.5" wide
- 18" seat
- Available in Red & Blue