



EM *obility*

Enhance Your Life
www.enhancemobility.com
800-296-4151

TRANSFORMER FOLDING SCOOTER



The Transformer Scooter is the original, revolutionary, electric folding scooter that opens and closes with a remote control. This popular travel scooter is so lightweight and compact, you can take it everywhere.

- Fully Assembled
- No Tools Required
- Fold & Unfold Quickly
- Fits in Any Vehicle
- Light Weight (54lbs)
- Easy to Carry
- Ideal for Travel
- 13.5 Miles Per Charge
- Armrests included
- Durable (300lb capacity)
- Adjustable Tiller



Available in 4 Colors




Mobie Plus

The most innovative folding scooter available. Its patented design allows you to easily fold and unfold in seconds. It's compact, lightweight and easy to use. The ideal scooter for travel, shopping and everyday use.

- Stylish Design
- Lightweight (54lbs)
- 300lb Weight Capacity
- Lithium Ion Battery
- Arm Rests Included
- Keyed Control Unit
- Range up to 20 miles*
- Ideal for Travel
- Folds/Unfolds in Seconds
- Available in Red & Blue
- Fits in Most Vehicles
- Adjustable Tiller
- Space Saving Design

*when traveling with spare battery

Compact, Lightweight, Folding – Perfect for Travel



TRIAXE SPORT



The unique TRIAXE Sport is a compact, portable, folding, electric scooter that can travel at speeds to 12mph and 35 miles on a single charge. It's the ultimate portable scooter that combines controlled speed and superb distance.

- Lightweight (44lbs)
- Aluminum Frame
- Rust Free
- Durable
- 350lb Weight Capacity
- Fast - Speeds to 12mph
- Up to 35 Mile Range
- Folds Flat for Transport & Storage
- Deluxe Basket Included

TRIAXE

SCOOTER



Available in 6 Colors

*Colors may vary



TRIAXE CRUZE

The TRIAXE Cruze offers the portability and convenience of a folding, travel scooter with the comfort and performance of a full size one. Why settle? when you can have it all!



- Lightweight
- Aluminum
- Rust Free
- Rear Wheel Drive
- Easy Paddle Acceleration
- Durable Aluminum
- 350lb Weight Capacity
- Electromagnetic Brakes
- Free Wheel Mode
- Speed Control
- Folds Flat in Seconds
- Great for Travel
- Basket Included
- Speeds to 6mph
- Up to 18 Mile Range



TRIAXE 
SCOOTER

HERCULES

LIFT


No More Struggling with A Heavy Scooter!

This ingenious lift folds compactly and makes lifting the Transformer, Mobie Plus and other similar scooters much easier. Simply unfold the lift and attach the power cord to the rear charging outlet on the scooter. The remarkable technology uses the battery power of the scooter to lift or lower the scooter in or out of the trunk of any vehicle. Easy to take along wherever you go.




- **Easy to Lift**
- **Motor DC 24V**
- **Lifts up to 66lbs**
- **Lightweight (29lbs)**
- **Lifts as High as 5 Feet**
- **Folds Compact & Portable**

Accessories Transformer & Mobie Plus



Accessories

TRIAXE Sport & TRIAXE Tour



This Page Left Blank for Notes



Enhance Mobility – the industry leader of compact, folding travel scooters and the force behind proven successes such as the Transformer Scooter, Mobie Plus, TRIAXE Sport and TRIAXE Tour. These scooters are designed to be light weight, portable and fun.

Our scooters are unique, not the typical mobility scooter. They appeal to a younger, more active audience. We believe an injury or disability shouldn't slow a person down. These compact, portable, travel scooters help to get people out and stay active. Our scooters take people further, faster and leave them with plenty of energy to enjoy their destination

All our scooters are approved for airline travel. Perfect for shopping, traveling, cruising and more. Your customers can take it with them wherever they go.

Our goal is to make the purchasing experience both easy and pleasant for our dealers. Shipping is always free (on purchases over \$89) within the continental USA.

Our dedicated service department and nationwide repair network are top-notch so you can rest assured knowing your customers will be taken care of after the sale.

To learn more about our dealer program, contact us at 800-296-4151 ext 704 or email sales@enhancemobility.com. Enhance Mobility – Enhance Your Life.



Please Note – all product specifications and colors are subject to change.