

**Its Freedom Baby!
YEAH!**

mojo
AUTO-FOLD



**REMOTE
OPERATED FOLDING**

**LIGHT
WEIGHT & COMPACT**

**13 MILE
BATTERY RANGE**

EMobility

Enhance Your Life



LED Display



Dual Lithium Ion Batteries



Detachable Transaxle

Model	MOJO Auto-Fold (MJ101)
Overall Length	39 3/4"
Overall Width	19"
Tires	Front Rear
	HoneyComb 6" Solid HoneyComb 7" Solid
Maximum Speed	Up to 3.75 MPH
Range	Up to 13 Miles
Turning Radius	55"
Weight	64 LB's
Color	Red, Aqua, Lemony Lime, Blue
Brake	Electromagnetic
Drive System	Rear Wheel Drive
Max Weight Cap.	Up To 300lbs
Seat Width	16.5"
Seat Height	23"
Battery	Dual(2) 24V 7.5AH Lithium
Motor	24V 120W
Battery Charger	29V 2AH
Folded Length	20.5"
Folded Height	27"
Heaviest Piece	Front Frame- 44 lbs



Enhance Your Life

For more Information...
www.EnhanceMobility.com
 1-800-296-4151
Sales@EnhanceMobility.com