



WHY SETTLE? YOU CAN HAVE IT ALL!

- **CONVENIENCE & PORTABILITY OF A FOLDING TRAVEL SCOOTER**
- **COMFORT & PERFORMANCE OF A FULL-SIZE SCOOTER**
- **DURABLE ALUMINUM CONSTRUCTION**
- **350 LB WEIGHT CAPACITY**
- **TOP SPEED OF 6 MPH**
- **LITHIUM ION BATTERIES- UP TO 18 MILES PER CHARGE**



Enhance Mobility

TRIAXE >>> CRUZE



Model	Triaxe Cruze (T3055)	
Overall Length	39"	
Overall Width	19.5"	
Tires	Front	8" Solid Tire
	Rear	9.5" Solid Tire
Maximum Speed	Up to 6 MPH	
Range	Up to 18 Miles	
Turning Radius	30"	
Weight	58 LB's	
Color	Fiery Red, Reflex Blue	
Brake	Electromagnetic	
Drive System	Rear Wheel Drive	
Max Weight Cap.	Up To 350lbs	
Seat Width	18" Standard	
Seat Height	23"	
Battery	24V 12AH Lithium Ion	
Motor	24V 180W, Brushed	
Battery Charger	24V 2AH	



Available Colors

Reflex Blue



Fiery Red



Enhance Your Life

For more Information...
www.EnhanceMobility.com
 1-800-296-4151
Sales@EnhanceMobility.com