

TRANSFORMER 2

FOLDING SCOOTER

Ultra Lightweight

Remote Operated Folding Scooter

Meets FAA Guidelines for Airline Travel



EMobility

Enhance Your Life



Specifications	Transformer S3026
Overall Length	38"
Overall Width	22" w/ Armrests
Tires	Front – 6x1.5" Rear – 7x2.36"
Maximum Speed	3.75mph
Range	Up to 13.5 miles
Turning Circle	55"
Weight (w/ battery)	55lbs
Color	Red, Blue, Yellow, Black
Brake System	Intelligent Regenerative Electromagnetic Braking
Drive System	Rear Wheel Drive
Ground Clearance	3 1/4 "
Weight Capacity	Up to 300lbs
Seat Width	17"
Seat Height (from ground)	22"
Battery	24V10Ah Lithium Ion
Motor	120W
Battery Charger	DC24V/2A
Folded Dimensions	18" x 22" x 25"



Enhance Your Life
 For more Information...
www.EnhanceMobility.com
 1-800-296-4151
Sales@EnhanceMobility.com