

**TRIAXE** >>>  
SCOOTER

# *Experience* the Triaxe Sport



**REDISCOVER ADVENTURE >**



Model	Triaxe Sport (T3045)	
Overall Length	34"	
Overall Width	20.5"	
Tires	Front	8" Solid
	Rear	9.5" Solid
Maximum Speed	Up to 12 MPH	
Range	Up to 35 Miles	
Turning Radius	33"	
Weight	46 LBs	
Colors	Red, Blue, Purple, Orange, Black, Silver	
Brake	Hand Controlled Disc Brakes	
Drive System	Front Wheel Drive	
Max Weight Cap.	Up To 350lbs	
Seat Width	13"Standard, 18" Optional	
Seat Height	20.5"	
Battery	48V 12AH, Or 48V 6AH	
Motor	48V - 500W	
Battery Charger	48V - 2AH	
Folded Length	34"	



**Available in Six Limited Edition Colors**



Enhance Your Life

For more Information...  
[www.EnhanceMobility.com](http://www.EnhanceMobility.com)  
 1-800-296-4151  
[Sales@EnhanceMobility.com](mailto:Sales@EnhanceMobility.com)